




# AseptoPeps



Powerful concentrated alkaline cleaner for milking equipment and milk tanks.

- ▲ Highly concentrated
- ▲ Free of chlorine and nitrogen
- ▲ "CompActive" powder

## APPLICATION

Pre-clean with water. Clean with **AseptoPeps**. Use 20 g of **AseptoPeps** powder per 10 litres of cleaning solution. Rinse with water of drinking quality. To be used in rotation with **HorolithPeps**. Do not prepare the solution in advance. Powder recipient must be clean and dry. Water temperature must not be below 35°C at the end of the cleaning.



## FEATURES AND BENEFITS

**AseptoPeps** is a strong alkaline product with boosted cleaning properties, which is suitable for all kind of milking equipment, including cooling tanks. It is able to eliminate all kind of milk deposits. The double concentrated and readily soluble powder leads to less packaging waste and transport.

The **Peps** products combine suitability to all water hardnesses with minimized application consumption of 0.2%.

## SPECIFICATIONS

AseptoPeps	
Physical form	White powder
Active substance	Caustic soda
pH at 20°C:	13.0 - 13.5
Density at 20°C (g/cm <sup>3</sup> ):	0.96 - 1.00
Storage temperature:	-20°C to +40°C

## LOGISTICS

Packaging	10 kg
Product code	2336410
EAN code	4028162336412

A Material Safety Data Sheet is available on request for professional users.  
A Product Data Sheet is available on request.

**Ecolab Europe GmbH**  
Richtstrasse 7  
CH-8304 Wallisellen  
www.ecolab.com