

IoShield® Spray

Double-performing sprayable shield to disinfect and protect teats after milking.



- ▲ High disinfectant efficacy of iodine
- ▲ Long lasting Shield protection
- ▲ Easy to apply and remove



APPLICATION

IoShield Spray® is a ready-to-use product. After each milking apply **IoShield Spray®** on the whole teat and do not wipe it. Keep the animals standing or walking for 5 min to ensure biocidal effect.

FEATURES AND BENEFITS

IoShield Spray® contains iodine, which due to its unique PVP carrier eliminates bacteria and disinfects the teats.

The quickly drying, polymerizing Shield of **IoShield Spray®** protects the teats between milking. It is used as part of a herd health plan for the prevention of mammary infections and provides protection from pathogens outside the udder and the environment.

IoShield Spray® provides prevention and protection irrespective of any cow, herd or environmental challenges. Well visible on the teats once applied, **IoShield Spray®** allows a quick and safe milking routine. Easy to remove, the formulation enables easier cleaning.



SPECIFICATIONS

IoShield Spray®	
Physical form	Dark brown liquid
Active substance	Iodine (PVP)
pH at 20°C	5.0 - 5.5
Density at 20°C	0.99 - 1.03
Storage temperature	+5°C to +25°C

LOGISTICS

Packaging	20 kg	58 kg	205 kg	1000 kg
Product code	2327310	2329140	2329150	2329160
EAN code	4028162327319	4028162329146	4028162329158	4028162329165

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com